## SDO<sup>®</sup> Original distributed by WELL & ABLE Sensory Dynamic Orthosis



Medigarments Ltd® Designed Around You


### PCP08 Sock Order Form

			Patient Deta	ails			
Order Detail	ed in order to process y	our order	Patient Reference No.: Is this their first SDO®?NoYes First Name:				
	-						
	Date Requi						
Order No.: Contact Name: Contact Phone No.: Email:							
			_ Date of Birth: Diagnosis:				
Hospital/Clinic: Delivery Address:			Attached with or	·	Consent form*		
				must be supplied when providi			
Post Code:			Clinician Number:				
Specification	IS						
Ankle Seam Angle	e: 45° (Standard)	90°					
Silicone Edging:	No	Yes					
Zip Position (	please specify)						
Legs:	Back	Front	Inside	Outside			
Themes							
1 Spaceman	2 Astronaut	3 Footballer	4 Striker	5 Goalie	6 Noir		
7 Silver Aztec	8 Pink Aztec	9 Leopard Flower	10 Lavender	11 Daisy	12 Ultraviolet		
13 Mermaid	14 Starlight	15 Marshmallow	16 Tribal	17 Blue Frenchie	18 Fawn Frenchie		
19 Cadet	20 Lance Corpora	I 21 Pilot Officer	22 Admiral	23 Major	24 General		
				g to the theme selected. e cannot proceed with ye			
Plain Coloure	ed Garments (B	ase fabric, panels, binding	g and thread are all the	same colour)			
White	Beige	Pink	Black	Navy Blue/Black*			
				*Navy base fabric, black panels, binding & thread			

When completed, please click: thrive@wellandable.com.au to email your electronic order form

# SDO<sup>®</sup> Original distributed by WELL & ABLE wearable therapy specialists

#### Sock Order Form



Important: All socks are supplied open toe unless closed toe is specified below. A closed toe option is not usually recommended.

#### **Circumference Measurements (C)**


16	С	Below knee at upper margin of tibial tuberosity		
17	С	Level with bulk of calf, 5mm less than 'at rest' circumference		
18	С	Level with mid-lower leg, at point of musculotendinous junction		
19	С	Level with upper margin of medial malleolus		
24	С	Foot and heel distal to the malleolus		
Instep	С	Around instep or waist of the foot		
26	С	Metatarsal heads		

#### Length Measurements (L)

21	L	Top of sock to upper margin of medial malleolus	
27	L	Upper margin of medial malleolus to sole of foot posteriorly	
28	L	From metatarsal heads to heel on medial border	
28A	L	From metatarsal heads to heel on lateral border	
29	L	From metatarsal heads to tip of toes medially. Required for closed toe sock only (not usually recommended)	
Zip	L	Required length of zip from top edge	

``.

#### Reinforcement Panels (Please specify)



© 2022 Medigarments Ltd. SDO® is the registered trademark of Medigarments Ltd.

Medigarments Ltd® Designed around you

Right (cm)

Right (cm)

Left (cm)

Left (cm)